


Reconstruction in CT and hints to the assignments

Jørgen Arendt Jensen

October 28, 2021


Center for Fast Ultrasound Imaging, Build 349
 Department of Health Technology
 Technical University of Denmark



Center for Fast Ultrasound Imaging
 Department of Health Technology


CT reconstruction repetition & hints

- Filtered backprojection algorithm
 - Radon transform
 - Filtered backprojection
 - Filters and their impulse responses
- Advise for the assignments
- Reading material: Prince & Links Chapter 6




From: W. A. Kalender; Computed Tomography, Publicis, 2005

Center for Fast Ultrasound Imaging, Department of Health Technology
 Technical University of Denmark



3


From: W. A. Kalender; Computed Tomography, Publicis, 2005



4

From: W. A. Kalender; Computer Tomography, Publicis, 2005

Parallel beam projection




Center for Fast Ultrasound Imaging, Department of Health Technology
 Technical University of Denmark

From: W. A. Kalender; Computed Tomography, Publicis, 2005

5

Parallel beam projection geometry



Center for Fast Ultrasound Imaging, Department of Health Technology
 Technical University of Denmark

6/x


6

Radon transform


$$p(x', \phi) = \int_{-\infty}^{+\infty} f(x' \cos \phi - y' \sin \phi, x' \sin \phi + y' \cos \phi) dy'$$

$f(x, y)$ – Attenuation image
 x', y' – Gantry coordinate system
 x, y – Patient coordinate system


7




8




9




10




11




12




13



14



15



16

Fourier Slice Theorem

$$P(\rho, \phi) = \int_{-\infty}^{+\infty} p(x', \phi) e^{-j2\pi\rho x'} dx'$$

$$= F(\rho \cos \phi, \rho \sin \phi)$$

$F(u, v)$ – 2D Fourier transform of image

ϕ – Gantry rotation

Filtered Back Projection (FBP)

$$\hat{f}(x, y) = \int_0^{\pi} \int_{-\infty}^{+\infty} |\rho| P(\rho, \phi) e^{j2\pi\rho x'} d\rho d\phi$$


$\hat{f}(x, y)$ – Reconstructed image

ϕ – Gantry rotation

x' – Detector position


Filtered backprojection

- Perform for all projections:
 - Make Fourier transform of projected data
 - Apply filter in Fourier domain
 - Make inverse transform
 - Back-project and sum with previous image



19

Influence from number of projections



20



Ram-Lak filter

Transfer function:

$$h(\rho) = \begin{cases} |\rho|, & |\rho| \leq B \\ 0 & \text{else} \end{cases}$$

Impulse response

$$h(k) = \begin{cases} B^2 & k = 0 \\ -\frac{B^2}{\left(\frac{\pi}{2}k\right)^2} & k \text{ odd} \\ 0 & k \text{ even} \end{cases}$$



Impulse response for Ram-Lak filter

- Impulse response of square wave:

$$h_{\text{square}}(x) = A2B \frac{\sin 2\pi B x}{2\pi B x}$$

- Ram-Lak impulse response: square wave – triangle:

$$h_{\text{ramlak}}(x) = 2BB \frac{\sin 2\pi B x}{2\pi B x} - B^2 \frac{\sin^2 \pi B x}{(\pi B x)^2}$$

Sampling of impulse response for Ram-Lak filter

- Sampling at $\frac{1}{\Delta x} = 2B, x = k \Delta x = \frac{k}{2B}$:

$$h_{ramlak}(k) = 2B^2 \frac{\sin 2\pi B \frac{k}{2B}}{2\pi B \frac{k}{2B}} - B^2 \frac{\sin^2 \pi B \frac{k}{2B}}{(\pi B \frac{k}{2B})^2}$$

$$\sin 2\pi B \frac{k}{2B} = \sin \pi k = 0$$

$$\sin^2 \pi B \frac{k}{2B} = \sin^2 \frac{\pi}{2} k = 1 \text{ (k odd) else } 0$$

- Gives:

$$h(0) = B^2$$

$$h(k) = -B^2 \frac{\sin^2 k \frac{\pi}{2}}{(k \frac{\pi}{2})^2} = -\frac{B^2}{(k \frac{\pi}{2})^2} \text{ (k odd)}$$

$$h(k) = 0 \text{ else}$$



Ram-Lak filter

Transfer function:


$$h(\rho) = \begin{cases} |\rho|, & |\rho| \leq B \\ 0 & \text{else} \end{cases}$$

Impulse response

$$h(k) = \begin{cases} B^2 & k = 0 \\ -\frac{B^2}{\left(\frac{\pi}{2} k\right)^2} & k \text{ odd} \\ 0 & k \text{ even} \end{cases}$$



Transfer function of filters – Ram-Lak



Center for Fast Ultrasound Imaging, Department of Health Technology
Technical University of Denmark

25/x

25


Shepp-Logan filter

Transfer function:

$$h(\rho) = \begin{cases} |\rho| \frac{\sin \frac{2\pi\rho}{4B}}{\frac{2\pi\rho}{4B}}, & |\rho| \leq B \\ 0 & \text{else} \end{cases}$$

Impulse response

$$h(k) = -\frac{8B^2}{\pi(4k^2 - 1)}$$




Center for Fast Ultrasound Imaging, Department of Health Technology
Technical University of Denmark

26/x

26

Shepp-Logan filter




Center for Fast Ultrasound Imaging, Department of Health Technology
 Technical University of Denmark

27/x

27

Hanning weighted filter




Center for Fast Ultrasound Imaging, Department of Health Technology
 Technical University of Denmark

28/x

28

Hanning weighted filter




Center for Fast Ultrasound Imaging, Department of Health Technology
 Technical University of Denmark

29/x

29


Filter transfer functions and impulse responses



30/x

30

Comparison between filters




Center for Fast Ultrasound Imaging, Department of Health Technology
Technical University of Denmark

31/x

31

Circular convolution




Center for Fast Ultrasound Imaging, Department of Health Technology
Technical University of Denmark

32/x

32

Circular convolution




Center for Fast Ultrasound Imaging, Department of Health Technology
Technical University of Denmark

33/x

33


Circular convolution – Shepp-Logan



Center for Fast Ultrasound Imaging, Department of Health Technology
Technical University of Denmark

34/x

34



35

Data for testing and validation

- Use data sets on web site
- Circular phantom for geometry test
- Shepp-Logan for orientation and quantitative data
- In-vivo images for Hounsfield units

TECHNICAL UNIVERSITY OF DENMARK

22485 Medical Imaging Systems

DTU Teaching Center for Fast Ultrasound Imaging Field II JAJ

CT data

On this page you can gain access to a database of CT images and projection data for the images, that can be used for making reconstruction algorithms. Data for some phantom objects and various programs for artificial phantoms and data projections are also given here.

The clinical images shown on these pages have been obtained from the Visible Human Project. They were created from scanning a human cadaver with CT and MRI scanners and then subsequently slicing the cadaver into 1 mm sections for taking photographs.

A further description of the data can be found at [here](#) and a description of the program can be found [here](#).

You can get to the different data pages by clicking on one of the images or the text below them.

Circle phantom Shepp-Logan phantom

Image of head Image of torso

Image of thorax Data for assignments

A slide showing Hounsfield units for a CT scan can be downloaded [here](#). The slide is taken from: E. Kwestel (ed.): Imaging systems for medical diagnostics, Siemens, 1990.

A table of Mass Attenuation Coefficients and Mass Energy-Absorption Coefficients can be found at the following reference:
<http://www.physics.nist.gov/PhysRefData/XrayMass/Coef/cover.html>

A radiologic image atlas can be found at:
<http://med.usfhs.mil/hadi/sona/homepage.html>

Biomedical Engineering Group DTU Teaching DTU Health Technology



Center for Fast Ultrasound Imaging, Department of Health Technology
Technical University of Denmark

36/x

36

Assignment data

- DICOM data from female patient
- All data available on the web
- Task is to find which slice it is



Data and program in:
 undervisning/k_22485_31545_billeder/ct_data/dicom_data

Center for fast ultrasound imaging, Department of Health Technology
 Technical University of Denmark

37/x

37

Reading DICOM data

```
% Set overall parameters
dir_name='DICOM/ST00001/SE00001/'; % Directory name
start_image=1; % First image in series
end_image=747; % Last image in series
frame_rate=50; % Frame rate for playing back the movie

% Set the dynamic range for the display
off_set=100; % Offset [Hu]
range=400; % Range to display [Hu]
map_values=128; % Number of gray level values
bone_off_set= -250; % Offset for showing the bones
bone_range=100; % Range for showing the bones

% Initialize figure
colormap(gray(map_values));
dicom_movie(end_image+1-start_image) = struct('cdata',[],'colormap',[]);

% Read information for the first images
file_name='IM00001';
info=dicominfo([dir_name, file_name]);
dx=info.PixelSpacing(1);
dy=info.PixelSpacing(2);
Y = dicomread(info);
[Nx,Ny]=size(Y);

% Make space for all the images
Y=zeros(Nx,Ny,end_image+1-start_image);
z_positions=zeros(end_image+1-start_image,1);
```

Center for Fast Ultrasound Imaging, Department of Health Technology
 Technical University of Denmark

38/x

38

```

% Loop through the images and read and display them

for i=start_image:end_image
    file_name=[ 'IM00', num2str(floor(i/100)), num2str(floor(rem(i,100)/10)), num2str(rem(i,10))];
    info=dicominfo([dir_name, file_name]);
    Y(:, :, i) = dicomread(info);
    image((1:Nx)-Nx/2)*dx, ((1:Ny)-Ny/2)*dy, (double(Y(:, :, i))+off_set)/range*map_values)
    xlabel('Distance [mm]')
    ylabel('Distance [mm]')
    pos=sprintf('%5.1f, %5.1f, %5.1f', info.ImagePositionPatient(1), ...
        info.ImagePositionPatient(2), info.ImagePositionPatient(3));
    z_positions(i)= info.ImagePositionPatient(3);
    text(-150, -150, ['Image ', num2str(i), ', Pos. (x,y,z) = ', pos, ' mm'], 'Color', [1 1 1])
    axis('image')
    drawnow
    dicom_movie(i)=getframe;
end

% Display the movie

movie(dicom_movie, 5, frame_rate);

```

Full script can be found at:

courses.healthtech.dtu.dk/22485/files/ct_data/dicom_data/display_dicom_images.m

on the page for the CT data: courses.healthtech.dtu.dk/22485/?ct_data/assign_data.html